

Grammage Checking Kit



GSM Round Cutter - (PRC-101)



Digital GSM Balance (SP-202)

GSM Round Cutter is a consistent tool useful for precise circular cutting of samples. Used to check the grammage of various sample. It cuts a sample with 100 sq. cm size rapidly and accurately.

Used in various Industries like paper and packaging, textile, films, foils etc.

- High Quality Aluminum pressure die casted main body
- Light weight and resilient
- Highly accurate test results
- Heavy Duty German Cutting Blades
- Available with Safety Latch Facility

Key Specifications:

Sample Area	100 sq.cm
Sample Diameter	112.8 mm
Handle material	Bakelite
Movement of Handle	¼ round to cut a sample

Highlights:

No. of Blades	4
Safety Latch	Yes
Material of Main Body	Die casted Aluminum
Sample size	100 sq.cm
Sample cutting pads	2 (Both side usable)

* Sample cutting thickness for paper & corrugated upto 4 mm.

Standard Accessories :

- Cutting pad pack of 2 nos. and Heavy Duty German Cutting Blades Pack of 4 nos.



GSM Balance Digital



Model No. SP-202

Table Top High Accuracy balance

- Capacity : 200gm
- L.C. : 0.01gm
- Pan size : 120 mm
- Display : high contrast LCD display
- Power : AC adaptor (Included) or 4 AA battery (not included)
- Calibration : 200 gm brass weight (Included)
- Linearity : ± 0.01 g
- GSM least count : 1 GSM = (weight x 100)

Optional : Air proof shield cover available on request.

* Other capacity 400gm, 600gm available on request

* Models with higher accuracy also available

Feature :

- Superior shock resistance & overload protection
- Software lock switch & integrated transportation lock
- Adjustable footrest with level indicator.

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-



Head Office : Presto Stantest Private Limited

I-42, DLF Industrial Area Phase-I, Delhi Mathura Road, Faridabad 121003, Haryana, India

P : +91 129 4272727 E : info@prestogroup.com

Branches :

Mumbai

503-504B, Fair link Centre, Opp. Fame Adlab
Behind Mongini's Off Link Road
Andheri (West), Mumbai 400058
P : +91 022 42957052, 42957053
E : west@prestogroup.com

Kolkata

93, 3rd Floor, Karnani Estate
209, A.J.C. Bose Road
Kolkata-700017. West Bengal
P : +91 033 40612012
E : east@prestogroup.com

Chennai

A-3, First Floor, Sri Sathya Sai Building
755, Poonamallee High Road, Kilpauk
Chennai- 600 010. Tamil Nadu
P : +91 044 43865260, 43838372
E : south@prestogroup.com

