



Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

- Faridabad • Kolkata • Mumbai • Pune • Ahmedabad
- Chennai • Bangalore • Hyderabad



Fume Hood



Model No. PFH-0125

Laboratory fume hoods are designed to protect laboratory personnel by preventing contaminants such as chemical vapors, dust, mists and fumes from escaping into the laboratory environment. Laboratory fume hoods also provide lab personnel with a physical barrier to chemicals and their reactions. Fume hoods are evaluated each year to verify their proper operation. They are generally set back against the walls to cover up the exhaust ductwork. Because of their shape, they are generally dim inside but have a provision for a light source. The front opening is made up of glass, able to move up and down on a counterbalance mechanism.

The Presto Fume Hood has heavy Stainless Steel Construction 304 Grade, a small S.S. Sink, Water Tap, Water outlet provided. A sliding glass door moves vertically with counter weight. Complete with fluorescent light. Centrifugal motor fan with S.S. fitted with the exhaust fumes will vented with PVC pipe outlet.



Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

- Faridabad • Kolkata • Mumbai • Pune • Ahmedabad
- Chennai • Bangalore • Hyderabad



Key Specification :

Test Chamber Size	2 Feet X 2feet
Blower Motor	1440rpm Single Phase
Main Door	Glass Viewing Single Hand Operated (Sliding Glass Door)
Work Table Material	Ceramic Glazed Tiles

Highlights :

Material	Out body Mild Steel powder coated & Inner body Stainless Steel.
Water Tap	Yes
Water Sink	Yes
Sample Viewing light	Yes
Finish	Powder coated Havel Gray & Blue combination finish and bright chrome / zinc plating for corrosion resistant finish

* Fume Exhaust point fitting & flexible pipe are not in our scope of supply.

* Water Inlet & water drain provision is not in our scope of supply.

Optional :

- Fume Hood with fully Stainless Steel body model is also available on request.

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-

Hot Air Oven Prima	Humidity Chamber	Cupping Machine	Salt Spray Chamber Korrox - III	Scratch Hardness Tester Motorised	Tensile Testing Machine

