

Wrap Reel Digital



Model No. PHW - 302

It is used to produce skeins of yarn of a pre determined length and number of turns for count and/or strength testing. It comes complete with yarn package stand and pretension device, fitted with Digital Preset Counter for auto stop after set cycles.

Special Magnetic brake motor stops at the exact count set on the counter to avoid extra round due to inertia.

- 5 position Reel.
- Circumference 1 meter
- Pretension device & yarn package stand.



Key Specification :-

Yarn Package Stand	5 Station
Wheel Circumference	1 Meter Standard (54 " also available on request)
Counter Type	Digital 7 Segment Display
Proximity Sensor	PNP
Counter Range	5 Digit (99999)
Counter Least Count	1 Rotation
Spread of each	20mm
Motor	100RPM \pm 10RPM ½ HP Break Motor
Power	Single Phase, 220V AC, 50HZ

Optional :

- Available in metric OR imperial system of units.

OUR OTHER PRODUCTS :-



Head Office :

Presto Stantest Private Limited

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India

P : +91 129 4272727 (100 Lines) E : info@prestogroup.com

• Faridabad • Kolkata • Mumbai • Pune • Ahmedabad • Chennai • Bangalore

www.prestogroup.com

