


Cobb Sizing Tester




Model No. PCS - 23

Cobb Sizing Tester is a consistent tool useful for checking the water absorptiveness (Cobb value) of a given test sample. Used to check the quantity of water absorbed by various sample like paper, cardboard, corrugated sample in a given time rapidly and accurately. PRESTO Cobb Sizing Tester is designed as per standards IS 1060, TAPPI –T441.


- High Quality rugged structure with corrosion resistant main body
- Light weight and resilient
- Fast Locking Mechanism
- Single handed hassle free operation
- Highly accurate test results
- High Grade stainless steel test area



Top View



Front View



Side View

Key Specifications:

Sample Diameter	130 mm (approx)
Depth of cylinder	60 mm
Diameter of rubber pad	145 mm
Water level indicator in cylinder	10 mm and 25 mm
Weight of roller	10Kg (approx)

Highlights:

No. of Roller	1
Locking screw	Yes
Material of Main Body	Stainless Steel / mild steel

Used in various Industries like paper and packaging, textile, films, foils etc.

Optional Features:

Combine with Digital Weighing balance for a complete Kit.

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-



Head Office : Presto Stantest Private Limited

I-42, DLF Industrial Area Phase-I, Delhi Mathura Road, Faridabad 121003, Haryana, India
 P : +91 129 4272727 E : info@prestogroup.com

Branches :

Mumbai

503-504B, Fair link Centre, Opp. Fame Adlab
 Behind Mongini's Off Link Road
 Andheri (West), Mumbai 400058
 P : +91 022 42957052, 42957053
 E : west@prestogroup.com

Kolkata

93, 3rd Floor, Karnani Estate
 209, A.J.C. Bose Road
 Kolkata-700017. West Bengal
 P : +91 033 40612012
 E : east@prestogroup.com

Chennai

A-3, First Floor, Sri Sathya Sai Building
 755, Poonamallee High Road, Kilpauk
 Chennai- 600 010. Tamil Nadu
 P : +91 044 43865260, 43838372
 E : south@prestogroup.com

