



Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

• Faridabad • Sonipat • Kolkata • Mumbai • Pune • Ahmedabad
• Chennai • Bangalore • Hyderabad

www.prestogroup.com



Muffle Furnace HMI - **PRIMA** SERIES



Test Report Graph



HMI Based
Touch Screen display

Muffle Furnace Prima Series is a unique test cabinet designed to detect the percentage of ash content in a particular test specimen like rubber, elastomers. It is Bench top model and is used to evaluate the amount of inorganic compound in a test material.

- Temperature sensing through J-type sensor.
- Highly accurate test results under uniform temperatures.
- PLC Control for Temperature.
- High Grade density imported Glass Wool Blanket insulation.
- Maximum Thermal Efficiency by insulation means.
- SSR based heaters.
- Real Time Data recording with temperature.
- USB option available for storing test data.
- Run time vs temperature graph.
- Over Temperature alarm system.





Presto Stantest Pvt. Ltd.

I-42, DLF Industrial Area Phase-1, Delhi Mathura Road, Faridabad 121003, Haryana, India
P : 9210 903 903, +91 129 4272727, 93111 24302 E : info@prestogroup.com

- Faridabad • Sonipat • Kolkata • Mumbai • Pune • Ahmedabad
- Chennai • Bangalore • Hyderabad

www.prestogroup.com



PRESTO Muffle Furnace comes with heaters and insulation maintaining heating efficiency and uniform temperature. It is applied in growing fields of Medical research, agricultural, industrial, metal treatment, ash determination & wool industries.

Key Specifications:

Size	<ul style="list-style-type: none">• 4x4x9 Inch• 6x6x12 Inch
Display	4.3 Inch HMI TFT screen display
Accuracy	± 5°C
Least Count/Resolution	1°C
Sensor	J-type/K-type as per selection of range
Power	20A, 220V AC, Single / Three phase, 50 Hz
Timer Range	Upto 999 hrs
Temperature Range	Ambient to 900°C / 1200°C as per requirement

* Selectable as per requirement. Other sizes available on request

Highlights:

Insulation	High grade ceramic insulation
Handle	Yes
Data Transfer	Through USB port
Material	Mild Steel (Outer Body) Ceramic Brick (Inner Body)
Finish	Powder coated Havel Grey & Blue combination finish and bright chrome / zinc plating for corrosion resistant finish

Optional Features:

Other Size and temperature range available.

Thank you customers for choosing us as your partners in growth !



*All the above logos are the sole property of their respective owners and are used purely for depiction purposes only.

OUR OTHER PRODUCTS :-

Top Load Tester	Wall Thickness Gauge	Gloss Meter Triangle	Eccentricity Tester	Moisture Analyzer MB 45	Preform Neck Cutter

